



Product Code	35264
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100405

## Description

OEM - Brand: Nihon-Kohden

N: 1

Sensor type: Adult Finger

Cable Length: 3.0 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Nihon Kohden: BSM-4101,4103 4111,4113 Life scope P BSM-5105, 5135 Life scope A BSM-2351, 2353 Life scope L BSM-2301, 2303 Life scope I, Life scope S