

SpO2 ADULT PROBE for NONIN - 0.9 m cable



Product Code	35266
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100092

Description

OEM - Brand: Nonin
N: 1
Sensor type: Adult Finger
Cable Length: 0.9 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: NONIN: 8500, 8600, 8700, 8800/9840; NONIN: 2500, 9600, 9700