

SpO2 NEONATAL "Y" PROBE for NOVAMETRIX- 3.0 m cable



Product Code	35268
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

EAN/UPC	6900234100580
---------	---------------

Description

OEM - Brand: Novametrix

N: 6

Sensor type: Neonatal Y-type

Cable Length: 3.0 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Novametrix: 505, 510, 511, 515B/C, 520, 7000, 7100; Oxypleth