

SpO2 ADULT PROBE for OHMEDA - 3.0 m cable



Product Code	35272
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

EAN/UPC	6900234100061
---------	---------------

Description

OEM - Brand: Ohmeda

N: 1

Sensor type: Adult Finger

Cable Length: 3.0 m

All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Datex-Ohmeda: 3740, 3710, 4700 OxiCap, 520RGM, Biox 3760, Ohmeda 3700, 3770, 3775, Rascal II, GE Medical System: 7000 series, Corometrics 511, TRAM A.R., Tram Modular x00SL Series, Mennen:Horizon XL800 (Ohmeda SpO2)