

SpO2 ADULT/NEONATAL WRAP PROBE for PHILIPS - 1.6 m cable



Product Code	35276
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100511

Description

OEM - Brand: Philips
N: 4
Sensor type: Adult/neonatal wrap
Cable Length: 1.6 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: M3, M4, MP20, MP30, MP40, MP50, MP70, MP80, VM4, VM6, VM8