

SpO2 ADULT PROBE for PHILIPS - 1.6 m cable



Product Code	35277
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100429

Description

OEM - Brand: Philips
N: 1
Sensor type: Adult Finger
Cable Length: 1.6 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: M3, M4, MP20, MP30, MP40, MP50, MP70, MP80, VM4, VM6, VM8