

SpO2 ADULT "Y" EAR PROBE for PHILIPS - 3.0 m cable



Product Code	35281
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100603

Description

OEM - Brand: Philips
N: 3
Sensor type: Adult Y-type ear
Cable Length: 3.0 m
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: M3, M4, MP20, MP30, MP40, MP50, MP70, MP80, VM4, VM6, VM8