

# SpO2 ADULT SOFT PROBE for PHILIPS - 3.0 m cable



Product Code	35283
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

CH-REP	No
--------	----

RDM (NSIS)	2621198
CND	C900301
EAN/UPC	6900234100474

## Description

OEM - Brand: Philips  
N: 2  
Sensor type: Adult soft-finger  
Cable Length: 3.0 m  
All products are supplied with manual in GB, FR, IT, ES.

Compatibility: Philips: 1205, M1205A, 78352C, 78834C, 78354C, M1020A, M1025A, M1025B, M1350B, M1772A/B, M1732A/B, M2475B