

EXTENSION CABLE for NELLCOR OXITECH 35262-3, 35265



Product Code	35305
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

CH-REP	No
RDM (NSIS)	83375
CND	Z12039080
EAN/UPC	6900234400062

Description

2.4 meter extension cable for Nellcor Oxitech SpO2 sensors.

It can be used with codes:

- 35262 Sensor type: Adult finger
- 35263 Sensor type: Adult soft finger
- 35265 Sensor type: Adult/neonatal wrap