

CABINET WITH ALARM FOR DEFIBRILLATORS - indoor use



Product Code	35331
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8034028012112

Description

METAL CABINET FOR DEFIBRILLATORS WITH ALARM - INDOOR USE - white

Door with opening at 180° with transparent acrylic window.
Suitable for wall mounting or use with pole.
Horn alarm at 100 dB

External size: 42.5x18xh 42.5 cm
Internal size: 40x15xh 40 cm