

# CABINET WITH THERMO AND ALARM FOR DEFIBRILLATORS - outdoor use



Product Code	35335
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8034028011573

## Description

METAL CABINET FOR DEFIBRILLATORS WITH THERMO REGULATION SYSTEM AND ALARM - OUTDOOR USE - yellow

Not suitable for marine or corrosive environments. Special coated cabinets resistant to mentioned environments available on request in 30 days.

Water tight welding on the sides, water and dust proof gaskets.

Made of carbon steel, thickness 15/10.

Adjustable internal shelf.

Door with transparent acrylic window.

With antitheft magnetic sensor that activates a 100 dB siren in case of non authorized withdrawal. 9V alkaline battery for alarm. Low charger indicator.

Power 230V for thermoregulation.

Suitable for wall mounting or use with pole.

External size: 43x21xh 48 cm

Internal size: 42x20xh 45 cm