

SP-10 POCKET SPIROMETER



Product Code	33536
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II A

RDM (NSIS)	2474966
CND	Z12150199
EAN/UPC	6945040100355
UMDNS	13674

Description

SP-10 POCKET SPIROMETER

Pocket, portable spirometer with large TFT screen to clearly display both numbers and flow charts.
Synchronously measure FVC, PEF, FEV1, FEV1%, FEF25, FEF25-75, FEF75.
Support diagnosis evidence by volume and flow rate chart display.
Built-in predicted values.
Free PHMS App to save data.
Built-in lithium rechargeable battery and battery power display.
Supplied with USB cable for PC connection, battery charger and software to calculate data and analyse trend.

Technical Specifications

Max volume: 10 l
Volume accuracy: $\pm 3\%$ or 50 ml whichever is greater
Flow range: 1 l/sec - 16 l/sec
Flow accuracy: $\pm 10\%$ or 0.3 l/sec whichever is greater
Parameters: FVC, FEV1, PEF, FEV1%, FEF25, FEF25-75, FEF75
CE 92/42 - EMC group I class B - protection type BF

Standard accessories

Internal lithium battery
USB cable
PC software (english language)
User manual: GB, FR, IT, ES, PT, SE, LT, SK, HU, RO