



Product Code	35360
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8056389239762

Description

METAL CABINETS FOR DEFIBRILLATORS - INDOOR/OUTDOOR USE

Metal cabinets for defibrillators - white

Alarm: no

Thermoregulator: no

Made of electro-galvanized carbon steel sheet, polyester powder coating high resistant to atmospheric agents.
Cabinets with rounded edges, hinged door with handle closure, accessory for plastic or metal seal and AED identification sticker.
Door with plexiglass glass, a rubber gasket on the ledge to ensure tightness and two sturdy painted hinges.
Holes for wall or totem application.

External size: 42x21xh 48 cm

Made in Italy.