

CABINET WITH THERMO AND ALARM FOR DEFIBRILLATORS - indoor/outdoor use



Product Code	35362
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	No
CH-REP	No
EAN/UPC	8056389239731

Description

METAL CABINETS FOR DEFIBRILLATORS - INDOOR/OUTDOOR USE

Metal cabinets for defibrillators - white

Alarm: yes

Thermoregulator: yes

Made of electro-galvanized carbon steel sheet, polyester powder coating high resistant to atmospheric agents. Cabinets with rounded edges, hinged door with handle closure, accessory for plastic or metal seal and AED identification sticker. Door with plexiglass glass, a rubber gasket on the ledge to ensure tightness and two sturdy painted hinges. Holes for wall or totem application.

Alarm

Cabinet case equipped with alarm and twilight internal lighting with 24-volt power supply.

Thermoregulation system and fan for air extraction

Extraction and recirculation fans are activated automatically when it exceeds 40 degrees C°, guaranteeing a temperature inside the cabinet in a range between 5 and 20 degrees C°. Exceeding temperature range is signaled by the flashing of red LEDs.

External size: 42x21xh 48 cm

Made in Italy.