

STAND for cabinets 35360-35362



Product Code	35363
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	No
--------	----

CH-REP	No
--------	----

EAN/UPC	8056389239724
---------	---------------

Description

METAL CABINETS FOR DEFIBRILLATORS - INDOOR/OUTDOOR USE

Stand for cabinets - white
Size: 42x21xh 100 cm

Complete with floor fixing holes and rear cover for cable passage

Made in Italy.