


BEAKER POSITIONING COVER 2x250ml for 35501-3



Product Code	35506
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279355062
---------	---------------

Description

Beaker positioning cover 2 x 250 ml for Branson 2800 ultrasonic cleaner (codes 35501-35502-35503).