

WIRE MESH BASKET for 35501-3



Product Code	35508
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279355086
PARAF	903115741

Description

Perforated tray for Branson 2800 ultrasonic cleaners.

- Compliant with:
- BRANSON 2800 M - with mechanical timer (code 35501)
 - BRANSON 2800 MH - with mechanical timer and heater (code 35502)
 - BRANSON 2800 CPXH - with digital timer and heater (code 35503)