


BEAKER POSITIONING COVER 2x600ml for 35501-3



Product Code	35509
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279355093
---------	---------------

PARAF	903115754
-------	-----------

Description

Beaker positioning cover 2 x 600 ml for Branson 2800 ultrasonic cleaner (codes 35501-35502-3)