


BEAKER POSITIONING COVER 3x250ml for 35510-2



Product Code	35516
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279355161
---------	---------------

Description

Beaker positioning cover 3 x 250 ml for Branson 3800 ultrasonic cleaner (codes 35510-35511-35512).