

BRANSON 5800 M ULTRASONIC CLEANER 9.5 l



Product Code	35520
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279355208
PARAF	903115842

Description

M Series (Mechanical timer): includes 60-minute mechanical timer and continuous operation.

Precision cleaning power
Ultrasonic cleaning is faster and safer than any other method.
It's powerful enough to remove heavy oils and waxes, penetrating enough to deflux circuit boards, yet gentle enough to clean tiny precision assemblies.

Ultrasonic sound waves moving through a cleaning solution create an effect called cavitation, the rapid formation and implosion of microscopic bubbles.

Heat enhances the process. The deep-cleaning action of ultrasonics removes even the most stubborn dirt and contaminations.

Branson: pioneer in ultrasonic cleaning
We have selected Branson because their cleaners incorporate years of experience in ultrasonic cleaning:

- easy to handle and safe to fill
- frequency higher than average for a deep cleaning
- less noise
- longer lifetime

All Branson cleaners are supplied with tank cover, and multilanguage manual (GB, FR, IT, ES, DE, PT, SE, NL, RU).

Recommended use with optional perforated tray.

High/low power control provides safer cleaning of delicate components. 100% power for normal or 70% for delicate applications.

Technical Specifications

- M series (mechanical)
- Timer: 60 minutes
 - Cover
 - Continuous mode
 - Sweep frequency
- BRANSON 5800
- Tank size: 290 x 240 x 150 mm
 - Overall size: 397 x 401 x 378 mm
 - Tank capacity: 9.5 l
 - Drain and valve kit: included
 - Weight: 6.4 kg
 - Frequency: 40 KHz
 - Power supply: 230-240 V