

# BEAKER POSITIONING COVER 4x600ml for 35520-2



Product Code	35529
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	802327935291
---------	--------------

PARAF	903115893
-------	-----------

## Description

Beaker positioning cover 4 x 600 ml for Branson ultrasonic cleaner 5510 (codes 33520-33521-33522)