

GLASS BEAKER 250 ml



Product Code	35540
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279355406
---------	---------------

PARAF	903115970
-------	-----------

Description

Glass beaker 250 ml for Branson ultrasonic cleaner:
- Branson 2800 (codes 35501-35502-35503)
- Branson 3800 (codes 35510-35511-35512)