



Product Code	35544
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279355444

Description

BRANSON PURPOSE is a cleaning solution developed for use with ultrasonic cleaners. Removes general soils, fingerprints, dust, light oils and greases.

Branson cleaning liquid has been specially developed to provide optimal transfer of acoustic energy.

Mainly used in: dentistry, medical and optical applications, laboratory, electronic, gold and silver work.

Technical Specifications

- Concentration: 10-12%
- Biodegradability: yes
- Phosphate free
- Concentration PH: 13.5
- Rinsability: good