



Product Code	35549
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279355499
---------	---------------

PARAF	903116022
-------	-----------

Description

This ultra-compact cleaner is very easy to use, efficient and economic. Ideal to clean small objects. Electronic, non adjustable timer with plastic cover and basket. No heating.

Internal size: 168 x 89 x 55 mm.

Supersonic wave 46,000 cycles per second.