

POSEIDONE BEDPAN WASHER



Product Code	35560
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	No
CH-REP	No

RDM (NSIS)	2303022
CND	Z12011301

Description

POSEIDONE BEDPAN WASHER

Designed to be used as a sink for organic substances and for washing and cleansing medical and sanitary devices such as: stainless steel or plastic urinals and bedpans, suction bottles, urine container.

Two washing programs: FAST cycle for urinals (max 130 seconds), NORMAL cycle for bed pans (max 300 seconds).

Front panel displays running status and the active alarm.

- electronic digital control
- light and acoustic warning
- function reset for changing the cycle
- function play/pause for the selected cycle
- instant stop of the cycle, in case of door opening
- made of stainless steel
- electronic check of electric pump

Multilingual manual: GB, FR, IT.

Made in Italy.

Technical Specifications

Automatic stand-by: at the end of the cycle
Maximum absorption: 2.8W in stand-by mode
Power supply - consumption: 230VAC - 50/60Hz - 400 W
Size: 536x735xh 930 mm