

TUTTNAUER COUNTER-TOP WASHER WITH DRYING SYSTEM



| | |
|---------------|----------------|
| Product Code | 35563 |
| Unit of sale | 1 pc |
| Minimum order | 1 |
| Type | Medical device |
| Class | II B |

| | |
|------------|---------------|
| UK-REP | Yes |
| RDM (NSIS) | 1872831 |
| EAN/UPC | 8056736494301 |
| UMDNS | 17671 |
| GMDN | 35424 |

Description

COUNTER-TOP WASHER with drying system

User-friendly counter-top washers to operate washing and disinfecting processes.

Washer chamber is made of high-quality stainless steel, AISI 316L.

Smart and efficient design, operated by a touch-keyboard and 3.5" LCD graphic display interface in 7 selectable languages (GB, FR, IT, ES, PT, DE, SE).

Three predefined operation programs and 37 customizable programs to choose from.

Both equipped with one dosing pump with flowmeters and level sensors.

Manual door with interlock and Integrated water softener.

Softener has the function of reducing the level of limestone contained in the incoming water.

- thorough drying using dedicated blowers
- integrated steam condenser
- real-time indication of water level and quality
- data transfer via USB and R232 with optional software
- two PT100 temperature sensors
- standard electrical connection (230V - 1 phase - 50Hz)

Supplied with standard washing basket, insert for vertical instruments, insert for 8 surgical cassettes, USB key with demo software* and manual in GB, IT (FR, ES, PT, DE, PL, RO, SE, RU on www.gimaitaly.com).

*Demo software has one week free trial, then the license must be purchased (code 35569)

Made in Italy.

Technical Specifications

Chamber capacity: 65 litres

Pump flow rate: 100 l/mn

Flow rate fan: 18 m³/h, 30 W

Power consumption: 2.9 kW

Internal tank: 0.6 l

Size - weight: 59.5x52x60 cm - 75 kg

Norms: EN ISO 15883-1/2; EN 1717

Operating voltage: 230 V - 50 Hz

Standard accessories

Cestello lavaggio standard (38,7x40,1x7,8 cm)

Inserto con 8 posizioni (18x37x10,5 cm)

Inserto per strumenti verticali (11,9x37x12,4 cm)