



Product Code	35570
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
--------	-----

RDM (NSIS)	2367975
CND	V9099
EAN/UPC	5214001621545

Description

T-CLEAN A TUTTNAUER ALKALINE DETERGENT - 5 litres

Highly concentrated and efficient alkaline detergent to clean surgical instruments, medical equipment, ophthalmological instruments, anesthetic and endoscopic equipment, laboratory glassware and similar applications. Formulated with alkaline agents, surfactants, highly stabilized enzymes, and powerful corrosion inhibitors, it makes the dissolution of organic residues such as fatty and protein deposits easier while effectively contributing to biofilm removal.

Recommended for heat-resistant and heat-sensitive instruments. It can be used in automatic washers, ultrasonic baths or immersion baths. It does not leave any residues. Used in soft or hard water. It does not foam, prevents corrosion and instrument discoloration.

Compatible with most materials: stainless steel, glass, ceramics, hard plastics, ebonite, etc.

Professional use.

Label: GB, FR, IT, DE. Medical device, class I.
MDR compliance (EU) 2017/745

Technical Specifications

Density: 1.02 g/cm³ at 20°C
pH: 10.00-10.50 at 20°C
pH (1% in deionized water): 10.20-10.80 at 20°C
Viscosity: <50 mPas at 20°C
Biodegradability: according to OCDE 301D

CHEMICAL COMPOSITION

Ingredients: < 5% non-ionic and anionic surfactants, chelating agent, corrosion inhibitor, pH regulator