



Product Code	35571
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2567268
CND	D0899
EAN/UPC	8009110030862

Description

GERMOCID ALKANET - tank 5l

Germocid Alkanet is a concentrated alkaline detergent which can effectively remove blood, protein, culture media, food residuals from medical instruments and medical devices. It is suitable for use in automatic washing equipment and for treatments in immersion or with ultrasound equipment. Its particular formula effectively removes dirt without altering the materials. It contains ecofriendly ingredients that prevent the sedimentation of organic residuals reducing the surface tension of the used solutions. Effective for the treatment of instruments in steel and stainless steel, optical equipment, conventional plastics and anesthesia instruments. For other types of materials (for example aluminum or light alloys) it is recommended to check the compatibility with the detergent. Use demineralized water for washing and rinsing.

Class I Medical Device.
MDR compliance (EU) 2017/745
Multilingual label: GB, FR, IT.
Multilingual instructions on www.gimaitaly.com: GB, FR, IT, ES, PT, DE, NL, IR, DK, SE, FI, LT, LV, EE, PL, CZ, SL, HU, SI, HR, RO, BG, GR, MT, RS, TR, RU, Arabic.

Made in Italy.

Technical Specifications

CHEMICAL COMPOSITION
100g of product contain: <5% amphoteric surfactants; phosphonates; enzymes; benzisothiazolinone