



Product Code	35575
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	2275515
CND	Z12011307
EAN/UPC	8059304700002

Description

TITANOX DRY STERILIZER - 3 I

Outer case made of stainless steel.
 Inner structure, grid, removable perforated tray and hinged lid completely made of stainless steel.
 Inner structure and lid with non-toxic silicon rubber gasket, insulated by high grade glass wool.
 Stainless steel armored heating element.
 Timer from 0 to 120 minutes, which disconnects the device at the end of the programmed time.
 Safety thermostat.
 Green pilot light (on the right) shows the line insertion.
 Green pilot light (on the left) points out the resistance operation.

Supplied with manual in GB, FR, IT (ES on request).

Made in Italy.

Technical Specifications

Capacity: 3 l
 Timer: from 0 to 120 minutes
 Probe thermometer: from 0° to 200°C
 Probe thermo-regulator: from 0° to 200°C
 Operating voltage: 230 V-50/60 Hz
 Overall dimensions: 400x295xh 170 mm
 Inside dimensions: 300x125xh 80 mm
 Classification: Class IIB, EU 2017/745 Medical Device