



Product Code	35576
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No
RDM (NSIS)	2275516
CND	Z12011307
EAN/UPC	8059304700019

## Description

### TITANOX DRY STERILIZER - 5 I

Outer case made of stainless steel.  
 Inner structure, grid, removable perforated tray and hinged lid completely made of stainless steel.  
 Inner structure and lid with non-toxic silicon rubber gasket, insulated by high grade glass wool.  
 Stainless steel armored heating element.  
 Timer from 0 to 120 minutes, which disconnects the device at the end of the programmed time.  
 Safety thermostat.  
 Green pilot light (on the right) shows the line insertion.  
 Green pilot light (on the left) points out the resistance operation.

Supplied with manual in GB, FR, IT (ES on request).

Made in Italy.

## Technical Specifications

Capacity: 5 l  
 Timer: from 0 to 120 minutes  
 Probe thermometer: from 0° to 200°C  
 Probe thermo-regulator: from 0° to 200°C  
 Operating voltage: 230 V-50/60 Hz  
 Overall dimensions: 500x325xh 170 mm  
 Inside dimensions: 400x160xh 80 mm  
 Classification: Class IIB, EU 2017/745 Medical Device