



Product Code	35581
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

RDM (NSIS)	2275524
CND	Z12011307
EAN/UPC	8059304700033
GMDN	35364

Description

TITANOX dry heat sterilizers

- STANDARD UNIT
- Capacity: 60 l
- Outside dimensions (mm): 705x450xh 475
- Inside dimensions (mm): 535x320xh 345

Inner structure and height adjustable inner shelves completely made of stainless steel.
Door with lock and key.

Standard unit: Double resistance for reaching high temperatures faster.

Supplied with manual in: GB, FR, IT, (ES on request) and 2 height adjust. trays

Technical Specifications

Timer: 0-120 minutes that switches off the oven at the end of the planned time

Probe thermometer: from 0° to 200°

Probe thermoregulator: from 0° to 200°

Operating voltage: 230 V - 50/60 Hz