



Product Code	35585
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	35046
CND	Z12011307
EAN/UPC	8023279355857

Description

TITANOX dry heat sterilizers

- THERMOVENTILATED UNIT
- Capacity: 120 l
- Outside dimensions (mm): 835x565xh 555
- Inside dimensions (mm): 670x420xh 415

Inner structure and height adjustable inner shelves completely made of stainless steel.
 Door with lock and key.

Thermovenilated unit: to guarantee the uniformity of the temperature in the heating room and to allow the circulation of the air, the device is provided with a ventilation motor that results always in operation, both during the cycle of heating and during the pauses of cooling.

- safety thermostat
- ventilator for forced ventilation
- double check up light bulb

Supplied with manual in: GB, FR, IT, (ES on request) and 2 height adjust. trays

Technical Specifications

Timer: 0-120 minutes that switches off the oven at the end of the planned time
 Probe thermometer: from 0° to 200°
 Probe thermoregulator: from 0° to 200°
 Operating voltage: 230 V - 50/60 Hz