



Product Code	35586
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	35041
CND	Z12011307
EAN/UPC	8023279355864

Description

TITANOX dry heat sterilizers

- DIGITAL UNIT
- Capacity: 20 l
- Outside dimensions (mm): 570x345xh 400
- Inside dimensions (mm): 405x255xh 210

Inner structure and height adjustable inner shelves completely made of stainless steel.
Door with lock and key.

Digital unit: adjustable aeration tunnels, placed on the bottom of the right side and on the top of the left side for hot air adjustment and circulation, which is useful for a quick cooling of the device.

- electronic running C.P.S.
- optimum heat distribution
- automatic indication for possible problems of feeding (in case of mains tension lack, the program will start again from the beginning)
- safety thermostat

Supplied with manual in: GB, FR, IT, (ES on request) and 2 height adjust. trays

Technical Specifications

Operating mode: continuous (stand by) or timed
Timer: from 0 to 250 minutes
Operating voltage: 230 V - 50/60 Hz