

RESISTOR 600 w for GIMETTE 21/28/50



Product Code	35639
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279356397

Description

Resistor 600W for sterilizer Gimette.

Specify the sterilizer model whit order:

- Code 35621 DRY STERILIZER GIMETTE 21
- Code 35628 DRY STERILIZER GIMETTE 28
- Code 35629 DRY STERILIZER GIMETTE 50