



Product Code	35641
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	1623433
CND	Z12011307
EAN/UPC	8054130053858

Description

MICRO STERIL QUARTZ BEADS STERILIZER 0.1 I.

Quartz beads sterilizers are ideal for sterilization of surgical tools, to eliminate the presence of bacteria, spores and micro-organisms. They ensure quick and effective sterilization at 240° C in 30 seconds.

Structure is entirely made of epoxy coated steel; the chamber is made of stainless steel.

Heat insulated with fiberglass.

Class II b.

Supplied with plug, pack of quartz beads (330 g) and manual in GB, FR, IT, ES, DE.

Made in Italy.

Technical Specifications

Operating voltage: 230 V
 Internal size: Ø 50xh 45 mm
 External size: 140x145xh 150 mm
 Capacity: 0.10 l
 Max. power: 70 W
 Timer: -