

# MICRO STERIL QUARTZ BEADS STERILIZER 0.2 l with timer



Product Code	35643
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	II B

UK-REP	No
CH-REP	No

RDM (NSIS)	1623436
CND	Z12011307
EAN/UPC	8054130053865

## Description

MICRO STERIL QUARTZ BEADS STERILIZER 0.2 l with timer.

Quartz beads sterilizers are ideal for sterilization of surgical tools, to eliminate the presence of bacteria, spores and micro-organisms. They ensure quick and effective sterilization at 240° C in 30 seconds.

Structure is entirely made of epoxy coated steel; the chamber is made of stainless steel.

Heat insulated with fiberglass.

Equipped with timer from 0 to 30 seconds.

Class II b.

Supplied with plug, pack of quartz beads (330 g) and manual in GB, FR, IT, ES, DE.

Made in Italy.

## Technical Specifications

Operating voltage: 230 V

Internal size: Ø 70xh 60 mm

External size: 145x190xh 155 mm

Capacity: 0.20 l

Max. power: 100 W

Timer: 0-30 seconds