



Product Code	35690
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279356908

## Description

### WATER DISTILLER

The autoclaves must run with distilled or demineralised water to avoid problems on electrovalves and on internal mechanism and white traces of limestone on the sterilized tools and devices.

Distiller produces 100% distilled water in a cheap and safe way.

It produces 4 litres in about 7 hours.

It is equipped with an automatic control thermostat with block system at the end of operation.

Supplied with 6 active coal filters, 1 detergent, 1 container for produced distilled water.

Operating Voltage: 220 V - 50-60 Hz.

Size Ø 20 cm x h 36 cm.

Weight: 3.4 kg.

It produces pure and clean water through a natural and effective process.

Water is brought gradually to 100 °C.

It eliminates:

- virus, bacteria, organic and inorganic chemical substances.
- heavy metals and contaminated agents settle in the bottom.
- impurities smells in the cooling coil.
- Water is further purified by the active coal filters present in the filter holder collector.

Multilingual manual: GB, FR, IT, ES, DE