

GASKET FOR H100 with serial number A,D,E,F,K,L



GUARNIZIONE / Gasket	
A D E F K L	} = 8 45 1045 JUN 0SB <input type="checkbox"/>
B C	} = 8 45 1045 099 0SB <input type="checkbox"/>
G	} = 8 45 1045 DE3 0SB <input type="checkbox"/>
H	} = 8 45 1045 MAG 0SB <input type="checkbox"/>
J	} = 8 45 1045 M99 0SB <input type="checkbox"/>

ETICHETTA MATRICOLA / Matriculation label

Model		
Code	1 11 1105 000		
Mtr.	99 ... A 0151	Year	1999
Max kPa	200	lt.	9
Max °C	130	Hz	50
	V 230~	A	F8
CE	W 1500		
Made in Italy			

Product Code	35747
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required

EAN/UPC	8023279357479
---------	---------------

Description

Gasket for Gima vertical autoclave H100 (cod. 35746-35656).

WARNING:

This gasket is only compatible for the autoclaves H100 with the following serial numbers:
A - D - E - F - K - L (see image, advise if you have different serial number).