


# COMPATIBLE PADS for defibrillator Nihon Kohden see also 55037



Product Code	33575
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
--------	-----

RDM (NSIS)	2474136
CND	C020401
GTIN	08054953346892
GMDN	45806

## Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market. Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Nihon Kohden
- Cardioline TEC 2202 - 2203 - 5500 - 7531 - 7600 - 7700 - JC-765V, TEC 5600 - 8300 with cable JC855V