

STERILIZATION BOX 1.5 l - plastic



Product Code	35770
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279357707
---------	---------------

PARAF	903119511
-------	-----------

Description

STERILIZATION BOX - 1.5 l (internal: 10x25.5xh 7.5 cm)

Useful plastic box, with removable internal perforated instrument tray, developed for cold sterilization of medical instruments.

Made in Italy.