


STERILIZATION BOX 3 I - plastic



Product Code	35772
--------------	-------

Unit of sale	1 pc
--------------	------

Minimum order	1
---------------	---

Type	No medical device
------	-------------------

Class	
-------	--

UK-REP	Not required
--------	--------------

CH-REP	Not required
--------	--------------

EAN/UPC	8023279357721
---------	---------------

PARAF	903119523
-------	-----------

Description

STERILIZATION BOX - 3 I (internal 18.3x29xh 9 cm) - autoclavable

Useful plastic box, with removable internal perforated instrument tray, developed for cold sterilization of medical instruments.

Made in Italy.