


COMPATIBLE PADS for defibrillator

Drager, Innomed, S&W, W-Allyn see also 55036



Product Code	33579
Unit of sale	kit of 2 pcs
Minimum order	1
Type	Medical device
Class	II B

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2468991
CND	C020401
EAN/UPC	8033003260524
GMDN	47055

Description

ADULT Compatible pads - disposable

A range of compatible pads BSI certified for patients >8 years (25 kg) for most defibrillators on the market.

Excellent adherence and adequate gel thickness.

Includes a couple of pads, cable and connector.

Maximum shelf life 30 months (24 months for 33584) from production date.

Latex free.

Compliance to regulations: Directive 93/42/EEC and amendments. Class IIb medical device.

ANSI/AAMI: DF-80; IEC/EN: 60601-1, 60601-2-4; ISO: 10993-1.

Multilingual envelope: GB, FR, IT, ES, PT, DE, NL, SE, GR, RU.

Made in Italy.

Compatible with:

- Artema Innomed
 - Cardio Aid 100 - 200 - 200B - 360B - XL - XL+
- Artema S&W
 - Cardio Aid LS - MC - MC+
- Drager
 - Cardiolog 2000
- Innomed
 - CardioAid 200B - 200 - 100 - 360B - XL+ - XL
- MRL
 - Lite - Lifequest
- Nihon Kohden
 - Cardiolife TEC 2202 - 2203 - 5500 - 7531 - 7600 - 7700 - 8251 all with cable JC-755V, TEC 5600 - 8300 with cable JC855V
- S&W
 - Cardio Aid LS - MC - MC+
- Welch-Allyn
 - AED 10 - AED20
 - PIC30 - PIC40 - PIC50