



Product Code	35806
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755779
CND	S01020101
EAN/UPC	8023279358063
GMDN	13735

Description

FLAT REELS for sterilization by steam or ethylene oxide gas

- FLAT REELS
- External: 200 m x 150 mm
- Internal: 200 m x 128 mm

Made of transparent, polyester-polypropylene film coupled to medical grade paper which permit steam and gas penetration during sterilization.

Steam and gas indicators provide clear colour changes when exposed to sterilization process.

Made in Italy.