

# ELINE GUSSETED REELS 100 m x 100 mm x 50mm



Product Code	35825
Unit of sale	box of 4 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755780
CND	S01020102
EAN/UPC	8023279358254
GMDN	13735

## Description

GUSSETED REELS for sterilization by steam or ETO - 100 m

- Gusseted reels
- External (m x mm x mm): 100x100x25
- Internal (m x mm x mm): 100x78x25

Made of coated film PET 12/PP 40 green and 60 g medical white paper.  
For sterilization of voluminous objects.

Made in Italy.