



Product Code	35827
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755780
CND	S01020102
EAN/UPC	8023279358278
GMDN	13735

## Description

GUSSETED REELS for sterilization by steam or ETO - 100 m

- Gusseted reels
- External (m x mm x mm): 100x200x30
- Internal (m x mm x mm): 100x178x30

Made of coated film PET 12/PP 40 green and 60 g medical white paper.  
For sterilization of voluminous objects.

Made in Italy.