

ELINE GUSSETED REELS 100m x 250mm x 65mm



Product Code	35828
Unit of sale	box of 2 pcs
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes

RDM (NSIS)	2755780
CND	S01020102
EAN/UPC	8023279358285
GMDN	13735

Description

GUSSETED REELS for sterilization by steam or ETO - 100 m

- Gusseted reels
- External (m x mm x mm): 100x250x32.5
- Internal (m x mm x mm): 100x228x32.5

Made of coated film PET 12/PP 40 green and 60 g medical white paper.
For sterilization of voluminous objects.

Made in Italy.