



Product Code	35830
Unit of sale	1 pc
Minimum order	1
Type	Medical device
Class	I

UK-REP	Yes
CH-REP	Yes
RDM (NSIS)	2755780
CND	S01020102
EAN/UPC	8023279358308
GMDN	13735

Description

GUSSETED REELS for sterilization by steam or ETO - 100 m

- Gusseted reels
- External (m x mm x mm): 100x400x40
- Internal (m x mm x mm): 100x378x40

Made of coated film PET 12/PP 40 green and 60 g medical white paper.
For sterilization of voluminous objects.

Made in Italy.