



Product Code	35846
Unit of sale	box of 12 pcs
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8001814292072

Description

STEAM TAPE INDICATORS

Tape indicator - 50 m: Steam Tape - autoclave
Application: Steam - moisture resistant
Width 19 mm

Gima range of tapes is specially formulated for packing purposes as well as to show visibly that the pack has undergone a sterilization cycle.

The backing of the different tapes is made of strong, semicreped, impregnated paper which gives excellent conformability during the application.

Sterile tapes unwind easily and can be torn by hand without difficulty.

Special adhesives formulated to withstand high temperature and humidity ensure excellent bond to paper, plastic, non-wovens, board, metal, glass and linen.

Each tape is individually packed with due date.

Steam tape:

for sterilization at 121°C - 10 minutes: from light to dark green

for sterilization at 134°C - 2 minutes: from light to medium / dark green

Multilingual instructions: GB, FR, IT, ES, PT, DE, NL, DK, FI, RO, HR, GR.

Made in Italy.