



Product Code	35853
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	3760243050039

Description

Bowie & Dick test ready for use.

- Stability of colour change.
- High sensibility, detect any anomaly
- Satisfactory: green uniform colour
- Not Satisfactory: Central red colour, Time and/or temperature and/or steam pressure uncorrect

The Bowie & Dick test pack is designed to control the proper functioning of steam autoclave and to test its ability to sterilise porous devices (compresses, fields, gowns, woven or non-woven, ...).

Operation:

The Bowie & Dick test pack Upstar is composed of a ready to use test pack including a test sheet placed in the middle of several cardboard layers. This test sheet is coated with a thermo-reactive chemical pigment changing color according to three parameters: temperature (134°C), time (3 min, 30 s), saturated steam.