



Product Code	35929
Unit of sale	1 pc
Minimum order	1
Type	No medical device
Class	

UK-REP	Not required
CH-REP	Not required
EAN/UPC	8023279359299

Description

GD-301 EVO SEALING MACHINE - complete

Manual heat sealer for medical pouches and reels according to DIN 58953-7:2008.

3 ceramic heater elements for a faster and more accurate sealing and improved infeed and outfeed guide to keep film flat (infeed) and prevent the film moving upward inside machine.

Available in 2 versions: compact without accessories (35928) or with sealer roll holder, cutter and plane (35929).

User friendly, reliable and extremely compact, ideal in dentistry, medical and veterinary laboratories.

Sealing bar 300 mm manually actuated by a lever.

Thermo controlled continuous heat sealing, to save consumables.

Compact vertical design to reduce space.

Integrated, removable roll-holder, can stand alone apart from the sealer.

Usable only with polypropylene film coupled to medical grade paper and/or plastic envelopes suitable for steam sterilization.

Technical Specifications

Flat sealing 12 mm
 Sealing length 300 mm
 Thermo controlled continuous heat sealing
 Built-in cutting device
 Fixed sealing temperature and time
 Manual actuation
 Power 230 V / 50-60 Hz, single phase
 Electrical input: 900 W
 Protection degree IP 20
 Sealer dimensions 452x248x146 mm
 Roll-holder dimensions 458x412x242 mm